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Sample collection for mycoplasma test

- Collect 1 ml of media to Eppendorf tube and store in -20°C until the test-PCR
- Anja Heikkinen (377B) takes care of two routine checks per year, otherwise, if needed, the researcher takes care of checking his/her own samples
- Ideally the supernatant should be collected 24 hours after changing the medium. If impossible, at least overnight. If the samples are collected after 48-72 hours of culture (without medium changes) you might get results that look like a "smear", because of DNA fragments from dead cells and debris.

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<https://notio.oulu.fi/fi/tilat-ja-palvelut>

Mykoplasmatestaus/Mycoplasma test:

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342B: Heta Merikallio, 348B, heta.merikallio@oulu.fi

357B: Miia Vierimaa, miia.vierimaa@oulu.fi

Instructions for working in cell culture

A dedicated lab coat must be used at any time when working in cell culture. Outside footwear should not be worn in the cell culture laboratory. Use shoes for only indoor use or shoe covers. Make a reservation for the hood beforehand. Reservation calendar works both in Outlook and O365 web calendar. Please contact Salla Kangas (salla.kangas@oulu.fi) to get access to the cell culture reservation calendar. If you are having problems with making reservations or access to the calendar, you can contact Salla.

Work aseptically:

- Gloves and a dedicated lab coat must be used at any time. Disinfect the gloves with ethanol prior to working.
- Wipe the microscope and other surfaces that are in contact with your cell culture dishes with ethanol before and after observing the cells.
- Laminar flow hoods are for sterile work. Keep them as empty as possible and do not overload them during working to ensure correct airflow to maintain sterility.
- Do not switch off the hood air flow.
- UV-lamps in the hoods are on for 1 hour between 5-6 AM.
- For routine cleaning of the hood, wipe the work surface with 70% ethanol before and after work, and during work especially after any spillage.
- Wipe the outside of the containers, flasks, plates, and dishes with 70% ethanol before placing them in the cell culture hood. To keep plastic ware and reagents sterile, they should be opened only inside the hood.
- Each hood has dedicated pipette controller and possibly micropipettes. Do not take pipettes out from their dedicated hood or cell culture room.
- Use only sterile filter tips in cell culture.
- Solid biological waste, (e.g. used pipette tips, cell culture dishes etc.), is collected into autoclavable bags. NOTE: Waste bin is not for liquids!
- Biological liquid waste is deactivated with 2% chlorine or virkon overnight. On the following day, pour the inactivated liquid waste down a sink and rinse with excess water.
- Glass pasteur pipettes have their own waste container. Also needles and other sharp waste has to be disposed in a specific container.
- Everyone is responsible of autoclaving their own glass and plastic ware. Everyone takes their own dishes out from cell culture.
- After using the vacuum aspirator, aspirate 2% chlorine/virkon and ethanol.
- Recharge the pipette controller when needed.

Cell culture waste:

- Solid biological waste, (e.g. used pipette tips, cell culture dishes etc.), is collected into autoclavable bags. NOTE: Waste bin is not for liquids!

- Biological liquid waste is deactivated with 2% chlorine or virkon overnight. On the following day, pour the inactivated liquid waste down a sink and rinse with excess water.
- Glass pasteur pipettes have their own waste container. Also needles and other sharp waste has to be disposed in their own specific container.
- After using the vacuum aspirator, aspirate 2% chlorine/virkon and ethanol. When the bag is full, replace it with new bag. Add chlorine/virkon into the used bag and on the following day, pour the inactivated liquid waste down a sink and rinse with excess water.

If you observe a contamination in your cell culture, dispose the contaminated culture and contact immediately cell culture responsible person Reetta Hinttala, phone: 050 428 6095, reetta.hinttala@oulu.fi

1. Dispose the contaminated culture and contact immediately cell culture responsible person Reetta Hinttala, phone: 050 428 6095, reetta.hinttala@oulu.fi
2. Document the contamination, if possible (e.g. take a picture of the dish/flask). After that, treat the dish with 2% chlorine/virkon.
3. Discard the potential sources of the contamination, such as culture media, PBS, etc. to avoid new contaminations.
4. Cell cultures that have been kept in same incubator as the contaminated dish are under constant observation. If no contamination is observed, the experiments can be finished as planned. Do not put new cultures in the incubator during that period. As soon as possible, the incubator is cleaned. DO NOT transfer the cells from this incubator to other incubators to avoid the risk to spread the contamination!

KEEPING THE CELL CULTURE LABS CLEAN

- UV-lamps in the hoods are on for 1 hour between 5-6 AM.
- When you use centrifuges, take care to leave it clean for the next user. Especially blood and other spills should be wiped away first wiping with water and then with ethanol.
- Testing for mycoplasma is done generally couple of times per year by Anja Heikkinen or by responsible researcher whenever needed.
- Biocidal usage according to own evaluation and needs in laminar flow hoods (please inform all other possible hood users).
- Patient-derived primary tissue cultures are concentrated in cell culture rooms 342B and 357B.
- Commercial and well characterized, clean cell lines are concentrated in an incubator on right side in room 341B, primary patient cells can be cultured in the incubator on the left hand side.
- Regular cleaning:
 - every 2 weeks
 - water baths
 - water baths inside the incubators
 - centrifuges
 - clean lab coat
 - every 4 weeks
 - clean lab coat

- wiping the laminars, also the lower bottom part

every 6 months

- incubators are sterilized